



ScaleBio™ Sample Fixation Kit

Laboratory Protocol

Scale Biosciences_Confidential_ScaleBio™ Sample Fixation Kit_RevA_09Mar23_Protocol

Table of Contents

<i>Introduction</i>	3
<i>Required Materials</i>	3
<i>Best Practices</i>	4
<i>Before You Begin</i>	5
<i>Buffer Preparation</i>	6
<i>Cell Preparation</i>	7
<i>Cell Fixation</i>	8

Introduction

The ScaleBio™ Sample Fixation Kit is intended for fixing cells prior to use in ScaleBio library preparation protocols. Fixation can be performed at different timepoints, as samples can be stored for up to 1 month at -80°C prior to use, reducing batch effects during library preparation.

Required Materials

Consumables and reagents contained within the ScaleBio Sample Fixation Kit (PN 2020001):

Kit Module	Consumable	Part Number	Storage Temp
ScaleBio Sample Fixation Kit – Module A (PN 2020002)	Wash Buffer	202100001	-20°C
ScaleBio Sample Fixation Kit – Module B (PN 2020003)	Fixation Regent	202110001	4°C

Consumables and reagents manufactured by other vendors:

Consumable or Reagent	Supplier	Part Number
1X PBS without calcium or magnesium	Various	Various
DEPC	Millipore Sigma	D5758-25ML
DMSO	Thermo	D12345
Methanol	Millipore Sigma	34860
Sterile, filtered, low retention tips for P1000, P200, P20, P10 pipettes	Various	Various
Sterile, filtered, wide bore tips for P1000, P200 pipettes	Various	Various
15-mL conical tubes	VWR	10025-286
50-mL conical tubes	Falcon	352070
1.5 mL LoBind Eppendorf tubes	Eppendorf	0030108418
Flowmi 40 µm filters	VWR	10032-802
Cell counting dye: AO/PI, Trypan Blue, YOYO-1, etc.	Various	Various

Required Equipment:

Item	Supplier	Part Number
Temperature-controlled centrifuges (for 15-mL & 1.5-mL tubes)	Various	
Vortex Mixer	Scientific Industries	SI-0236
P1000 pipette	Various	
P200 pipette	Various	
P20 pipette	Various	
P10 pipette	Various	
Cell counting instrument: Hemocytometer, Nexcelom Cellometer K2, etc.	Various	
Cell counting slides	Various	

Best Practices

For general laboratory practices:

- Calibrate and service pipettes every 12 months to ensure accurate sample volume transfer at each step.
- Store all reagents at the storage conditions recommended by the supplier.
- Unless otherwise specified, thaw reagents on ice.
- Unless otherwise specified, vortex reagents.
- Open Fixation Reagent packing in a chemical fume hood.
- Handle Fixation Reagent, DEPC, and Methanol in a chemical fume hood.
- Never reuse pipette tips or tubes.
- Use wide-bore tips for pipetting cell mixtures.
- Keep pipette tip boxes, reagent containers, and sample tubes closed when not in use.
- Wear suitable protective clothing, eyewear, and gloves.

For RNase-free sample processing:

- Use low-retention, RNase-free pipette tips and low-binding reaction tubes to prevent adsorption to plastic surfaces.
- Routinely wipe work surfaces with RNase AWAY to remove RNases, and with a 10% bleach cleaning solution to remove DNA amplicon contaminants.
- Wear disposable gloves and change them frequently.
- Keep pipette tip boxes, reagent containers, and sample tubes closed when not in use.
- Routinely wipe work surfaces with a 10% bleach solution.

Before You Begin

- Bring tubes of Fixation Reagent to room temperature.
- Fully chill 100% Methanol on ice.



Note: Centrifugation of 15-mL conical tubes and 1.5-mL microcentrifuge tubes are performed at 4°C. Bring centrifuges that accommodate these two tube formats to 4°C.

- Place 1X PBS on ice.
- Thaw Wash Buffer at room temperature, then invert Wash Buffer to ensure it is fully mixed. Place the Wash Buffer on ice.
- Determine the speed setting for the vortex mixer that will be used during the fixation process, using a 15-mL conical tube containing 500 μ L of water. The speed of the vortex mixer should be set such that the height of the 500 μ L liquid reaches the 5 mL mark on a 15-mL conical tube as shown below:

Figure 1. Setting the vortex speed

500 μ L water without vortexing



500 μ L water with vortexing



Buffer Preparation



Important: Prepare buffers using a vortex mixer in a chemical fume hood.

Materials:

- Fixation Reagent (1 tube per sample)
- DMSO (50 μ L per sample)
- Ice-cold 100% Methanol (2 mL per sample)
- 15-mL conical tube (1 tube per sample)

Procedure:

1. Prepare Complete Cell Fixation Solution
 - a. **In a chemical fume hood**, reconstitute Fixation Reagent by adding 50 μ L DMSO to each tube. Fixation Reagent is lyophilized at the bottom of the tube and appears as a white pellet.



- b. Dissolve tube contents by repeatedly vortexing at high speed with intermittent brief spins to bring contents to the bottom of the tube. Repeat this process until all solids are fully dissolved; this may take up to several minutes. Ensure that all solids are fully dissolved before proceeding.
- c. Briefly spin down the Fixation Reagent tube to bring contents to the bottom of the tube.
- d. Prepare Complete Cell Fixation solution by combining ice-cold 100% methanol with reconstituted Fixation Reagent in a conical tube according to the volumes provided in the table below:

Reagent	1 Sample	2 Samples	4 Samples
100% Methanol	2 mL	4 mL	8 mL
Reconstituted Fixation Reagent	50 μ L	100 μ L	200 μ L

- e. Vortex for 10 seconds to mix and place on ice. Use within 6 hours.

Cell Preparation



Important: We recommend use of wide-bore pipette tips, and gentle pipette-mixing, for preparing cell suspensions to maintain sample quality. A total of 2.5 million or fewer cells can be processed per tube of Fixation Reagent. Cells must be washed out of their preparation medium with 1x PBS without calcium or magnesium and resuspended in a maximum volume of 500 μ L prior to fixation.



Note: If starting samples are $\leq 400,000$ cells, please contact your FAS for suggestions on protocol modifications for lower cell numbers.

Materials:

- Ice-cold 1X PBS without calcium or magnesium
- Sterile, filtered, wide-bore pipette tips
- 15-mL conical tubes
- 40 μ m Flowmi Cell Strainers
- Cell counting dye (Trypan Blue, YOYO-1, AO/PI, etc.) and required instrumentation (hemocytometer, Countess, Nexcelom Cellometer K2, etc.)
- Counting slides

Procedure:

1. Obtain cells from culture or prepare from frozen.
2. Wash cells with ice-cold 1X PBS.
3. Centrifuge at 500 x g for 5 min at 4°C.
4. Remove supernatant taking care to not aspirate the cell pellet.
5. Resuspend cells with 500 μ L 1X PBS.
6. Determine cell concentration using a viability dye and a hemocytometer, Nexcelom Cellometer K2, or similar cell counting equipment. For accurate cell counting, use ≥ 2 μ L of cell suspensions and appropriate dilution factors recommended for your cell counting method.
7. If cell suspensions appear clumpy, strain cells to a new 1.5-mL tube through 40 μ m Flowmi cell strainer and repeat step 6 to recount the cell suspension.
8. Transfer 2.5 million or fewer cells to a 15-mL conical tube and place on ice.
9. Bring the volume of the cell suspension up to 500 μ L with 1X PBS.

Cell Fixation



Duration: 45 minutes

Materials:

- Flowmi strained and counted single cell suspension in a 15-mL conical tube
- DEPC (20 μ L per sample)
- Complete Cell Fixation Solution (2 mL per sample)
- Wash Buffer (one bottle per sample)
- Sterile, filtered, wide-bore P1000 pipette tips



Note: All steps below are done on ice with centrifugations done at 4°C.

Procedure:

1. Immediately before use, add 20 μ L of DEPC per 2 mL of the Complete Cell Fixation Solution (prepared on page 6, step d). Briefly vortex to mix.
2. Using the settings determined in Figure 1, vortex the cells while adding 2 mL of the Complete Cell Fixation Solution + DEPC dropwise to the tube.



Caution: Adding the fixative too quickly to the cells can result in cell clumping/incomplete fixation.

3. Incubate the tube on ice for 15 minutes.
4. Vortex the cells while adding 5 mL of Wash Buffer dropwise to the tube.
5. Centrifuge the tube at 500 x g for 5 min at 4°C.
6. Carefully remove supernatant without disturbing the pellet, leaving ~50 μ L of residual volume as shown:



7. Gently flick the cell pellet until the cells are fully resuspended in the residual volume in the tube.
8. Using a wide-bore pipette tip, add 1 mL of Wash Buffer while rinsing down the sides of the tube. Transfer the fixed cells to a new 1.5-mL tube.
9. Centrifuge the tube at 500 x g for 5 min at 4°C.
10. Carefully remove supernatant without disturbing the pellet, leaving ~50 µL of residual volume.
11. Gently flick the cell pellet until the cells are fully resuspended in the residual volume in the tube.
12. Add 100 µL of Wash Buffer and gently flick the tube until cells are fully resuspended.
13. Determine cell concentration using a viability dye and a hemocytometer, Nexcelom Cellometer K2, or similar cell counting equipment. For accurate cell counting, use ≥ 2 µL of cell suspensions and appropriate dilution factors recommended for your cell counting method.



Note: Fixed cells may settle at the bottom of the tube. To ensure even distribution of cells, flick the tube 10-15 times until pellet has dispersed before counting cell suspensions.

14. If fixed cells will not be used within the same day as preparation, freeze cells at -80°C.



Frozen cells can be stored at -80°C for up to 1 month prior to use. We suggest storing aliquots at a concentration of between 5,000-25,000 cells/µL.

The ScaleBio Single Cell RNA Sequencing Kit is protected by U.S. Patent Nos. 10,626,442; 10,982,256; 11,512,341; 11,634,752; 11,566,278